



## Quick Reference Chart TPLO Clover indications



Ref.	V1861 (L) V1862 (R)	V2063 (L) V2064 (R)	V2061 (L) V2062 (R)	V3061 (L) V3062 (R)	V3063 (L) V3064 (R)	V3065 (L) V3066 (R)
Description	<b>TPLO Clover Micro</b> Length: 26.6 mm Thickness: 1.5 mm	<b>TPLO Clover Mini XS</b> Length: 37 mm Thickness: 2.2 mm	<b>TPLO Clover Mini</b> Length: 44 mm Thickness: 2.5 mm	<b>TPLO Clover Medium</b> Length: 55.3 mm Thickness: 2.8 mm	<b>TPLO Clover Large</b> Length: 60.5 mm Thickness: 2.8 mm	<b>TPLO Clover Broad</b> Length: 79 mm Thickness: 3.0 mm
Weight Kg.	3-5	up to 10	10-15	up to 35	up to 50	50-70
Blade	12	12-15	18	21-24	24-27	27-30

The values in this table are provided only as a guide. The implant choice has to be based on surgeons experience. Surgeons must also consider patients attitude, temperament and containment conditions.



**Intrauma S.p.A.** - [www.intrauma.com](http://www.intrauma.com) - Tel. +39.011.953.94.96 - [orders@intrauma.com](mailto:orders@intrauma.com)  
**Intrauma America LLC** - P. (201) 420.0110 - [orders@intrauma.com](mailto:orders@intrauma.com)

